









Realistic Fluorophore Separation Capacity Analysis

Accounting for Spectral Overlap and Lifetime Distributions



REALISTIC SEPARATION CAPACITY ANALYSIS

Instrument Specs:

- Wavelength resolution: 10 nm
- Lifetime resolution: 0.5 ns
- Polarization states: 2

Theoretical vs Realistic:

- Theoretical bins: 30 (λ) × 9 (τ) × 2 (pol) = 540
- Effective spectral bins: 28 (not 30!)
- Effective lifetime bins: 6 (not 9!)
- Realistic capacity: 8 fluorophores

Key Limitations:

- Spectral overlap: Emission spectra are broad (FWHM = 35-65 nm, not single wavelength)
- Lifetime distributions: Not single values (σ = 0.1-0.4 ns variation)
- Crosstalk: Overlapping signals require unmixing

Practical Recommendations:

- Conservative panel: 3-5 fluorophores
- Moderate panel: 5-8 fluorophores
- Aggressive panel: 8-8 fluorophores (requires advanced unmixing)

Best Strategy:

- Maximize wavelength separation (most powerful)
- Use lifetime for same-wavelength species (e.g., NADH free/bound)
- Polarization adds modest 2× benefit

Example Well-Separated Panel (5 fluorophores):

- DAPI (461 nm, 2.5 ns)
- GFP (509 nm, 2.5 ns) - 48 nm separation
- YFP (527 nm, 3.0 ns) - 18 nm + lifetime
- mCherry (610 nm, 1.5 ns) - 83 nm separation
- Cy5 (670 nm, 1.0 ns) - 60 nm separation